M93UB Landmark puncture pad I Ultrasound puncture pad I

CVC Insertion Simulator II

Caution | Don't mark on the model and other components with pen or leave printed materials contacted on their surface.

Ink marks on the models will be irremovable.

Instruction Manual



Table of contents

Before you start Manufacturer's note P.1 Set includes P.2 Dos and Dont's P.3
 [I . Landmark puncture pad II] Parts names and functions ····· P.5 Preparation ···· P.6-P11 After a session ···· P.12-P13
[II. Ultrasound puncture pads II] Introductory ultrasound training block "REAL VESSELS" Parts names and functions P.16 Preparation P.17-P19 After a session P.20-P22 Storage and maintenance
[III. Transparent cannulation block] Parts names and functions · P24 Preparation · P.25-P26

Before you start

Manufacturer's note

Manufacturer's note

The simulator is designed for training in CVC procedures.

Any other use, or any use not in accordance with the enclosed instructions, is strongly discouraged. Kyoto Kagaku Co., Ltd. cannot be held responsible for any accident or damage resulting from such use. Please use this model carefully and refrain from subjecting it to any unnecessary stress or wear. Should you have any questions or concerns regarding use of this simulator, please contact the distributor you purchased from or Kyoto Kagaku Co., Ltd.

Features

The simulator comes with 3 kinds of training pads for relevant area and an introductory ultrasound training block.

■ Landmark puncture pad II

Three possible accesses: subclavian, internal jugular and supraclavicular.

Carotid artery pulsation is palpable.

A guide wire can be inserted to SVC.

Success and failure confirmation

Jugular/artery puncture can be known by the color of the fluid you collect.

Pneumothorax: your syringe will collect air only.

Guide wire (catheter) placement and Mislodging can be confirmed by finding the wire in the

Ultrasound puncture pad II

Both internal jugular and axillary veins are accessible under ultrasound scanning. Anatomically correct vein and artery relationship. Veins collapse under light pressure of the probe.

Transparent cannulation block

Facilitates three-dimensional anatomical understanding. As an effective training tool for developing guide wire insertion skills.

Introductory ultrasound training block REAL VESSELS

An introductory training block to acquire basics of ultrasound-guided puncture.

Before you start

Set includes

Set includes

Before your first use, ensure that you have all components listed below.



Α.	male upper torso manikin		
Β.	transparent anatomical block1		
C.	skin for cannulation training 1		
D.	landmark puncture pad II · · · · · · · · · 1		
E.	. ultrasound puncture pad II · · · · · · · · · 1		
F.	 introductory ultrasound training block ··· 1 "REAL VESSELS" 		
G.	irrigation bottle · · · · · · · · · · · · 1		

Н.	plastic jar · · · · · · 2
Ι.	Cotton Swab for Coloring each 1 pack (red / blue)
J.	air bulb · · · · · · 1
K.	50ml syringe · · · · · · · 1
L.	sample syringe · · · · · · · · · · · · · 1
M.	sample needle · · · · · · · · · · · · · · · · 1
	Instruction Manual

Before you start

DOs and DON'Ts

▲ DOs and DON'Ts

DOs

Use new needles for training. Recommended needle size: 23G or thinner

Durable transparent cannulation pad is highly recommended for insertion training in order to maintain the puncture pad longer, although it is possible to demonstrate the full procedure with puncture pad.

Handle with care.

The materials for the Simulator are a special composition of soft resin. Handle them with the utmost care at all times.

Clean the surface of the pad with dry or wet soft cloth and put talcum powder before storage. Clean the upper torso manikin by dry cloth.

Store the training set at room temperature, away from heat, moisture and direct sunlight.

DON'Ts

Never use the supplied needle for anything other than the simulator.

Do not peel the skin from the puncture pads.

Do not push the tube tip at the bottom of the landmark puncture pad.

Never wipe the simulator with thinner or other organic solvent.

Don't mark on the simulator and components with pen or leave any printed materials in contact with their surface. Ink marks on them will not be removable.

Do not store the simulator and pads with any water left inside.

Do not use broken or bent needles for training.



When you use the puncture pads; after aspirating do not return the fluid from syringe into simulator. It may damage the vessel tube and may lead to leakage.

Landmark puncture pad ${\rm I\!I}$



Parts names and functions	P.5
Preparation	
1. Attach the air bulb · · · · · · · · · · · · · · · · · · ·	P.6
2. Set the pad on the torso	P.6
3. Preparation of red fluid ·····	P.7
4. Fill the artery tube with red fluid · · · · · · P. 7-	P.8
5. Preparation of blue fluid	P.9
6.Fill the vein tube with blue fluid ••••••••••••••••••••••••••••••••••••	v .10
After a session	
 Discharge fluid from the pad from the artery tube with red fluid 	P.11
2. Discharge fluid from the pad from the vein tube with blue fluid ••••••••••••••••••••••••••••••••••••	v.12
3. Detach the puncture pad from the simulator F	? .13

Find the how-to-use video on our website

http://www.kyotokagaku.com

Parts names and functions



Attach the air bulb Place the pad into the torso

Attach the air bulb





Carotid artery pulsation can be made by grabbing the bulb.

Attach the pulsation unit to the connector tube at the rear wall of the landmark pad. Ensure to hold the connector tube by your finger tips so that the tube is not pushed toward the pad wall.

Attach the air bulb

▲ Caution

Do not put the pad on flat surface or push the connector for arterial pulsation unit, as this may damage the pad.

Set the pad on the torso



2





Putting the pulsation unit through the opening in the right shoulder of the torso, place the pad in the cavity. Put the tubes through the hole of the head side. Open the cut on the left shoulder and put the tubes into the cut.



Make sure that the tubes are not folded or tucked between walls of the pad and the body torso. The simulator doesn't work properly when the tubes are folded.



3

4

Preparation of red fluid Fill the artery tube with red fluid

Preparation of red fluid



Pour 200ml water into the plastic jar.Put the simulated blood(swab type/red) into the jar and stir the water sufficiently to prepare the simulated blood.

▲ Caution

Take care not to drop simulated blood on clothes as blood stains can be very difficult to remove. Don't drink the simulated blood.

Fill the artery tube with red fluid



1. Put the torso upright.

A Caution



2. Open the Clamp.



Before start filling the fluid. Ensure that the clamp is open. Otherwise exsessive pressure may cause damage to the pad.



Fill the artery tube with red fluid



3. Fill the 50ml syringe with red fluid and connect the syringe tip to the connector at the lower end of artery tube (upper-back tube from the right shoulder).



4. Then push the piston slowly until the red fluid reaches to the height of the clamp.

Insert and turn clockwise







8

5. Close the clamp.

6

Preparation of blue fluid Fill the vein tube with blue fluid

Preparation of blue fluid



Pour 200ml water into the plastic jar.Put the simulated blood(swab type/blue) into the jar and stir the water sufficiently to prepare the simulated blood.

▲ Caution

Take care not to drop simulated blood on clothes as blood stains can be very difficult to remove. Don't drink the simulated blood.

Fill the vein tube with blue fluid



Put the torso upright.

A Caution



Open the clamp.



Before start filling the fluid. Ensure that the clamp is open, otherwise exsessive pressure may cause damage to the pad.



Fill the vein tube with blue fluid



Syringe

Fill the vein tube with blue fluid

Right subclavian artery tube

To pad



3. Fill the 50ml syringe with blue fluid and connect the syringe tip to the connector at the lower end of vein tube (lower-front tube from the right shoulder).

Insert and turn clockwise

Detach



4. Then push the piston slowly until the blue fluid reaches to the height of the clamp.
Ensure that all tubes are fully filled with the fluid and no air bubbles remain inside.
By tilting the torso carefully, air can be removed.



5. Close the clamp.





1

Discharge fluid from the artery tube

Discharge fluid from the artery tube



 Place the torso upright.
 Connect the empty 50ml syringe to the connector at the lower end of the artery tube.



2. Open the Clamp of the artery tube.







 Pull the syringe slowly and drain the fluid.
 To clean the tube, fill the tube with uncolored water and discharge it again.

Discharge fluid from the vein tube





1. Connect the empty 50ml syringe to the connector at the lower end of the vein tube.



2. Open the clamp of the vein tube.







 Pull the syringe slowly and drain the fluid. To clean the tube, fill the tube with uncolored water and discharge it again.

Detach the puncture pad from the simulator

3 Detach the puncture pad from the simulator



Open the cut of left shoulder and remove the tubes. Pull out the puncture pad from inferior side.



Remove the puncture pad from the torso.



Disconnect the pulsation unit from the pad. The tube can be disconnected by pushing the stopper ring toward the connector.





Be careful NOT to push the connector tube tip toward the puncture pad. The pressure may cause damage to the soft structure of the puncture pad.

Ultrasound puncture pads ${\rm I\!I}$



Introductory ultrasound training block "REAL VESSELS"	
1. General information	P.15
2. Preparation	P.15
3. Training and after a session · · · · · · · · · · · · · · · · · · ·	P.15
Ultrasound puncture pad II	
Parts names and functions	P.16
Preparation	
1. Place the pad into the torso •••••••••••	P.17
2. Fill the artery tube with water •••••••• P.17-	-P.18
3. Fill the vein tube with water ••••••••••••••••••••••••••••••••••••	-P.19
After a session	
1. Discharge water from the artery tube •••••••	P.20
2. Discharge water from the vein tube	P.21
3. Detach the puncture pad from ••••••••••••••••••••••••••••••••••••	P.22
4. Maintenance of the skin of the pad •••••••••••	P.22

Find the how-to-use video on our website

http://www.kyotokagaku.com

REAL VESSELS

Introductory ultrasound training block "REAL VESSELS"

General information

This block facilitates training in basics of ultrasound guided punctures, before moving onto trainings with the anatomical type ultrasound pad.

Features

2 simulated vessel lines: straight and curving. Both lines are embedded with slope to represent from shallow vessel to deep one.

Vessel wall yields under pressure of a needle tip.

Training items

How to get correct ultrasound images.

Probe manipulation. Basics for ultrasound-guided vessel access.

Preparation



Peel the protective sheet carefully.





Pour water into the wider slit in the container and fill it up to the line on the wall. Ensure that shallower end of vessels are fully under the water surface.

Training and after a session



Spread some ultrasound gel on the surface of the block and start training. When the needle tip is in the vessel, your syringe can collect water.

Add water to the container slit as necessary.



After use, wash off the gel completely with running water and dry well. Then, replace the protective sheet before storage.



3

Do not leave paper, cloth or other materials except for the included protective sheet contact with the block surface. Such items may stick to the surface and lead to damage the block.

Parts names and functions



1

Place the pad into the torso Fill the artery tube with water

Place the pad into the torso



1. Put the tubes through the hole of the head side.



2. Open the cut of left shoulder and put through the tubes. Place the pad into the torso.



Make sure that the tubes are not folded or tucked between walls of the pad and the body torso. The simulator doesn't work properly when the tubes are folded.



(17)

2

3

Fill the artery tube with water Fill the vein tube with water

First, fill the artery tube with water



 Connect the water-filled 50ml syringe to the connector at the lower end of the artery tube (upper and rear side tube at the shoulder).
 Screw in the tip of the syringe to the connector. Advance the piston slowly until the water surface reaches height of the tube clamp.





 When the tube has been filled by water, close the tube clamp.
 Take off the syringe from the tube.



Then, Fill the vein tube with water



1. Put some water (around 5cm height) in the irrigation bottle.

Connect the tip of the tube from the bottle to the head end of the vein tube.



 Connect the syringe filled with water to the connector at the lower end of the vein tube.
 Screw in the tip of the syringe clockwise to the connector until it locks.

3

Fill the vein tube

Then, Fill the vein tube with water



3. Push the piston of the syringe slowly until the water flows out into the bottle. Take off the syringe from the tube. Keep the bottle connected while the simulator



4. Lay the simulator and start your training session.



When the vessel tube is empty, ultrasound image will not be shown properly.

19

Discharge water from the pad

Discharge water from the artery tube



1. Put the torso upright and connect the empty 50ml syringe to the lower end of the artery tube.





2. Open the tube clamp.





3. Drain water by pulling the piston of the syringe. Detach the syringe from the tube and discharge the collected water.



Ensure that all the water is discharged from the pad after use. Do not store the pad while any water remains inside as it may lead to malfunction or damage the pad.

2

Discharge water from the pad

Discharge water from the vein tube



 Put the torso upright and disconnect the tip of the tube from the irrigation bottle.
 Connect the empty 50ml syringe to the lower end of the vein tube.
 Drain water by pulling with the syringe.



2.Detach the syringe and discharge the water.





3.Tilt the body torso toward the syringe to aspirate remaining water in the tubings.



Ensure that all the water is discharged from the pad after use. Do not store the pad while any water remains inside as it may lead to malfunction or damage the pad.



3

Detach the puncture pad from the simulator Maintenance of the skin of the pad

Detach the puncture pad from the simulator



1. Open the cut of left shoulder and remove the tubes. Pull out the puncture pad from inferior side.



2. Remove the puncture pad from the torso.



Do not lift the pad by pulling the vessel tubes as this may cause breakage in the embedded tube walls.

4 Maintenance of the skin of the pad



1. To remove the ultrasound gel on the skin, wipe the skin surface gently by soft cloth, or soft wet close. If stickiness persists, wipe with soft cloth with small amount of alcohol for disinfection.



 The skin of the ultrasound pad is not fixed and may come off if you hold the pad by skin only. In such case, replace the skin not to leave air bubbles under it. Air bubble or layer between skin and soft tissue will impair the quality of scanning image.



Never use thinner or organic solvent for cleaning. It may cause breakage in the simulator skin. Do not touch the surface of the skinless pad, or expose the surface to air for a long time. Dust attached to the surface may weaken its adhesion and lead to deterioration in scanning images.

Transparent cannulation block



Pa	arts names and functions	P.24
G	eneral information	P.25
Ρ	reparation	
	Remove the pad currently set on the torso ······	P.25
	Set up the cannulation block	P.26
	Detach the pad from the torso	P.27

Find the how-to-use video on our website

http://www.kyotokagaku.com

Parts names and functions



General information Remove the pad currently set on the torso

General information

The transparent cannulation block is an effective educational model to facilitate understanding of the relevant anatomical structures and is applicable to trainings in cannulation by interchanging with the puncture pads.

Training items

A Caution

- Anatomical understanding
- Learn the appropriate depth and angle of the needle for each approach.
- Handling and manipulation of the catheter (guide wire).
- Simulate the steps of the procedures.

- Do not make any puncture at any site besides the prepared openings.
- The transparent block is not designed to be filled with fluid.
- Do not pour any fluid or water to the openings on the cannulation block.
- Please handle the skin sheet with utmost care. Excessive strain may cause breakage.



Set up the cannulation block

2 Set up the cannulation block



Place the transparent cannulation block to the cavity.



Now you can practice cannulation by observing the catheter through the transparent block.

Cover the pad with the skin sheet



Cover the cannulation block with the skin sheet. Put the skin sheet over the pad so that the slits of the skin sheet fit with the openings on the pad.



Now, the simulator is ready for cannulation training. The slits in the skin sheet allow trainees to remove it with the catheter inserted in the pad.



Do not make any puncture at any site besides the prepared openings. The transparent block is not designed to be filled with fluid. Do not pour any fluid or water to the openings on the cannulation block. Please handle the skin sheet with utmost care. Excessive strain may cause breakage.

Detach the pad from the torso





Remove all needles and catheters from the simulator. Remove the skin.



Remove the pad holding its frame.

* If you are moving onto next training session using other kind of training pad, set up the new pad following the instruction on :

Set up the landmark puncture pad	P.6
Set up the ultrasound puncture pad	P.17

Trouble shooting

L U us

Quick check-up before calling the customer service. Use the table if you have problems using the simulator. Look in this section for a description of the problem to find a possible solution. (TEL + 81-75-605-2510)

Trouble	Possible Reason	What to Do
l cannot fill/discharge	The body torso with the pad is laid down.	Put the torso upright.
the fluid to/from the vessel tubes properly.	One or some vessel tubes from the pad are folded.	Straighten the tubes.
	The puncture pad is worn out.	Order new pads for replacement.
	The puncture pad is worn out.	Order new pads for replacement.
Heavy leakage from puncture area.	Pool of fluid is formed in the pad.	Straighten the tubes. Withdraw all fluid from the tubes and squeeze the remaining fluid from the soft tissue. Dry the pad naturally.

L landmark

Carotid pulsation	The pulsation unit (air bulb and tube) is not connected to the pad.	Connect the pulsation unit to the connector at the bottom of the landmark pad.
does not work.	The tube of the pulsation unit is folded in the pit.	Make the tube straight.

Uus		
	The vessel tubes are not filled by water/fluid. Or, air bubbles are formed in the vessels.	Fully fill the vessels with water.
Ultrasound image is unclear.	Air layer is formed under the skin sheet.	Reattach the skin sheet removing air bubbled and layers.
	The puncture pad is worn out. Order new pads for replacement.	Order new pads for replacement.



Ink from this manual can stain the model - do not let it touch the model and its components.

For inquiries and service, please contact your distributor or KYOTO KAGAKU CO., LTD.



■ Main Office and Factory (World Wide)

KYOTO KAGAKU co., ltd

WEB•www.kyotogagaku.com E-MAIL•rw-kyoto@kyotokagaku.co.jp

15 kitanekoya-cho Fushimi-ku Kyoto 612-8388, Japan TEL: +81-75-605-2510 FAX: +81-75-605-2519 North and South American regions:



WEB•www.kkamerica-inc.com E-MAIL•info@kkamrica-inc.com

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA TEL: +1-310-325-8860 (Toll-free in North America: 877-648-8195) FAX:+1-310-325-8867

The contents of the instruction manual are subject to change without prior notice.

No part of this instruction manual may be reproduced or transmitted in any form without permission from the manufacturer. Please contact manufacturer for extra copies of this manual which may contain important updates and revisions. Please contact manufacturer with any discrepancies, typos, or mistakes in this manual or product feedback. Your cooperation is greatly appreciated.

Instruction manual is in the inside pocket of the carrying case. PLEASE **READ** THE **MANUAL** CAREFULLY **BEFORE USE**.

