Caution:

Don't mark on the phantom with pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable.

PH-2 Whole Body Phantom PBU-50

Instruction manual

Contents

Please read before training		
General information · · · · · · · · · · · · · · · · · · ·	P.1~	P. 3
Preparation		
Assembly manual ·····	P.4~	P. 9



Whole body phantom PBU-50

Please read

The PBU-50 is life-size human phantom with a life-size synthetic skeleton embedded in radiological soft tissue substitute.

There are no metal parts that would impair realistic imaging. The joints, knees, elbows, shoulders and hip-joints, are flexible and able to maintain the necessary positions.

The phantom is separable into 10 parts facilitating a variety of application. A positioning stand for the head part comes with the set.

Specifications

Lengt	ded application: Positioning phantom for plain radiography h: approx.165cm nt: approx. 50kg
Can b	e dismantled at: shoulders elbows, hip joints, knees, neck
Embe	dded organs: lungs with vessels, heart, liver and kidneys
	ne product has two 3mm dia holes made in the chest wall to prevent the damage caused v atmospheric pressure change that may occur in transit.
Mova	ble joints and their range of movement
9	houlders: rotate 360 degrees to back and forward, 180 degrees side-ways
(Cubital joints (bilateral): bend inward up to 90 degree
(Coxae (bilateral): rotate forward up to 90 degrees, then rotate outward up to 45 degrees each.
	Articulatio genus (bilateral): bend up to 90 degrees.
Mate	ials: Human tissue substitute (urethane, epoxy, etc.,)
Acces	sories: Positioning head stand 1 piece
	Flat head screwdriver 1 piece
	Sample images 1 set
	Hand positioning belt 1 piece
Repla	cement parts: screws and connection tapes for Patella.
I	umbers of each organ: iver 70
	idney 30 source and the second se

Please read

1 Please note

- 1) Bone trabeculae are observed only at hands.
- 2) Contours of the bones are slightly stranger than actual human x-ray image, to facilitate clear understanding.
- 3) Basically the phantom is designed to be radiographed under same setting as human body, though, the joints areas may require some adjustment to have better image.

▲ Caution

- Do not make the phantom to take the position unintended by the manufacturer. It may cause the breakage in the phantom if the joints are forced to bend or rotate to the directions or degrees of angles outside of designed movable range.
- 2) Do not remove the plastic protectors on the phantom's hands. The protector prevents the phantom fingers from breaking.
- 3) Do not attempt to lift or carry the phantom if the head or limbs are attached to the trunk. Detach the joints before moving to avoid the excessive force on the joints.
- 4) When the training with the phantom is not made for a while, pay attention to storage condition (posture) of the phantom so as to avoid the burden to the connetion area. For example, condition the head has been resting on the pillow for a while can cause the breakage to the neck connection area.
- 5) Do not fasten the screws too tight.
- 6) The phantom is heavy and could cause damage not only to itself but to people or other objects. When you need to handle the assembled phantom, always do it by two or more people with utmost care.
- 7) Don't mark on the phantom with pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable.

2

PBU-50 Components

Parts #	Parts name	QTY
1	Head	1
2	Neck-head connection plate	1
3	Trunk	1
4	Right upper arm	1
5	Left upper arm	1
6	Right forearm	1
7	Left forearm	1
8	Right thigh	1
9	Left thigh	1
10	Right lower leg	1
11	Left lower leg	1
12	Right patella	1
13	Left patella	1
14	Connector for shoulder A (Right or left)	2
15	Connector for shoulder B (Right or left)	2
16	Connector for right elbow (lateral side)	1
17	Connector for right elbow (medial side)	1
18	Connector for left elbow (lateral side)	1
19	Connector for left elbow (medial side)	1
20	Connector for right hip joint	1

Parts #	Parts name	QTY
21	Connector for left hip joint	1
22	Supporting plate for right hip joint	1
23	Supporting plate for left hip joint	1
24	Connector for right knee	1
25	Connector for left knee	1
26	Screw bolt for neck connection M12	2
27	Screw bolt for neck connection M10×15	2
28	Screw bolt for shoulder connection Screw bolt for parts # 22.23 M 8×10	10
29	Screw bolt for elbow joint M 6×60	4
30	Screw nut for elbow joint M 6	4
31	Screw bolt for hip joint M8×35	6
32	Screw bolt for knee joint M12×20	4
33	Connection tape for the patellae	2
34	Height adjustment panel	1
35	Base panel	1
36	Neck connector	1
37	Neck connector fixture	1
38	Screw for fixture	1
39	Screw bolt M6×18	2
40	Hand positioning belt	1



4

* The circled numbers in the assembly images are parts numbers on the components table on page 3.

< NECK AND HEAD >

1. Remove the two bolts for the neck connection (part #26) and remove the neck head connection plate (part #2) from the trunk.



Screw bolt connection M12×20



[Neck sub-assembly]

2. Screw the neck head connection plate (part #2) on the head (part #1).







Screw bolt for neck connection M10×15×2

(27)

[Connection of the head and the trunk]



3. Mount the assembled parts (parts #1 and #2) to the trunk and secure each side with a screw bolt (part #26).



< ELBOWS >



forearm

7

Screw bolt for 29 elbow joint M6×60×4 Screw nut for 30





< SHOULDERS >



1. Remove the two clear bolts (part #28) and remove the external shoulder connector parts A (part #14) and connector parts B (part #15).





(lateral side)

2. Then attach the upper arms (parts #4 and #5) at the ball joint to the trunk and place over the two connectors for the shoulder (parts #14 and #15). Secure the connectors with the clear bolts (part #28). Follow these steps for both left and right side of the trunk.







14)

(15)





.......... Screw the bolts after confirming the position of the screw hole. Be careful not to mix up part #4 and 5. The shape of right and left arms are different.

6

shoulder

<knees>

1. Connect the thigh (#8, 9) and lower leg (#10, 11) at the knee joint. Insert the dowel of the upper thigh through the hole of knee joint of lower thigh. Set the connector for knee (#24, 25) to the medial side then screw the knee joint bolt (#32) to fix it.



2. After fixing the knee joint, bend the knee and insert the connection tape of the patella (#12, 13) to the guide aperture of the thigh.



7





<HEAD AND STAND >

1. Mount the head phantom to the neck connector and fix it with two screws.





27) Screw bolt for neck connection M10×15×2



Then, insert the end of the neck connector (#36) in the slit in the panel and fix it with the fixture (#37).









- 3. The phantom is ready.
 - Adjust the height of the neck as needed, by loosen the fixture and slide it through the slit. Tighten the fixture again at the designated height.





Don't mark on the phantom with pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable.

Main Office and Factory (World Wide)

KYOTO KAGAKU co., ltd

WEB•www.kyotogagaku.com E-MAIL•rw-kyoto@kyotokagaku.co.jp

15 kitanekoya-cho Fushimi-ku Kyoto 612-8388, Japan TEL: +81-75-605-2510 FAX: + 81-75-605-2519 **North and South American regions:**



WEB•www.kkamerica-inc.com E-MAIL•info@kkamrica-inc.com

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA TEL: +1-310-325-8860 (Toll-free in North America: 877-648-8195) FAX:+1-310-325-8867

2020.06