US-17 Rheumatism Hand Phantom

Caution:

Don't mark on the phantom with a pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable.



Thank you for choosing our Rheumatism Hand Phantom.

This phantom is designed for a training of ultrasound imaging that is helpful for early diagnosis of articular rheumatism. For effective training, use ultrasound devices that is under appropriate control.

Features

- Training of ultrasound imaging of bones of finger joints and articular rheumatism symptoms.
- Hyperplasia of synovial membrane and condition of increase of synovial fluid are reproduced.

Before Using

Components and handling precautions

Including items

Check the following including items before using.

- A) Rheumatism Hand Phantom1
- B) Carrying Case 1
- C) Talcum Powder 1
- D) Instruction Manual.....1







Cautions	
This phantom is made of soft material, so	Immediately after training, wipe out gel on the
strong impact caused by fall so on may cause	phantom with wet tissue etc.
the damage.	*Never use any organic solvent such as
	thinner to wipe out gel.
Don't mark on the phantom with a pen or leave	The color of the surface may change when used
printed materials contacted on its surface.	for a long time. This color change does not affect
Ink marks on the phantom will be irremovable.	ultrasound image.
To avoid the damage of the phantom, please	After use, keep the phantom away from heat,
handle the phantom with care.	humidity and direct sunlight. These may
Do not bend, pull or screw the fingers of	cause damage of the phantom.
phantom.	
Do not pick the surface layer of the	
phantom with your fingers.	

Training

Preparation, Training and Cleaning

1) Preparation

This phantom is designed to allow training in scanning fingers from both the palm side and the back side of the hand.

To avoid the damage of the phantom, please lay a soft rug under the phantom while the training.



Apply gel for ultrasound to the probe.

2) Training

Scanning the phantom with a transducer.

Training points

- Setting for ultrasonography of rheumatism and holding of probe.
- Positioning
- Selection and operation of a transducer.
- Visualization of synovial membrane hyperplasia.
- Interpretation and classification of ultrasound images.



3 Cleaning

Wipe out gel on the phantom.

*Dried residue of gel may affect ultrasound images and cause damage of the surface of the phantom.





Don't mark on the phantom with a pen or leave printed materials contacted on its surface. Ink marks on the phantom will be irremovable



Ultrasound image of a joint of the phantom



KYOTO KAGAKU co., ltd

Kyoto Kagaku Co., Ltd Head Office.

15 Kitanekoya-cho Fushimi-ku Kyoto, Japan 612-8388 Telephone:+81-75-605-2510 Fax:+81-75-605-2519 Email:rw-kyoto@kyotokagaku.co.jp



Kyoto Kagaku America Inc. 3109 Lomita Boulevard, Torrance, CA 90505-5108,USA Telephone:1-310-325-8860 Fax:1-310-325-8867 Email:rw-kyoto@kyotokagaku.co.jp

Kyoto Kagaku Co., Ltd. European Representative Office

De-Saint-Exupery-Str.10, 60549 Frankfurt, Germany Telephone:069-5060-28160 Fax:069-5060-28170

2019.09