Ultrasound Compatible Lumbar Puncture/ Epidural Simulator

Lumbar puncture training with ultrasound guidance

Anatomically correct puncture block offers excellent ultrasound images







FEATURES

INAGES

- 1 | Durable and replaceable parts for repetitive practice
- 2 | Adjustable to an upright or lateral position.
- 3 | Translucent block makes the needle trace visible.
- 4 | Interchangeable parts for other landmark and fluoroscopic procedure training

SKILLS

| Ultrasound-guided lumbar puncture

- | Ultrasound-guided epidural anesthesia
- | CSF collection and CSF pressure measurement



SPECIFICATIONS Size: W33 x D21 x H30 cm W13 x D8.3 x H11.8 inch

SET INCLUDES

- 1 lumbar region model
- ultrasound lumbar puncture/epidural block
- 1 ultrasound lumbar region skin cover
- 2 lumbar region support bases
- (upright position / lateral position)

1 irrigation bag syringe

- 1 stand for irrigation bag
- instruction manual

DESCRIPTIONS

MATERIALS Soft resin, polyurethane elastomer Latex free

REPLACEMENT PARTS

11348-190 ultrasound lumbar puncture /epidural block 11348-230 ultrasound lumbar region skin cover for M43E

RECOMMENDED DEVICES Ultrasound machine and needles are not included.

KYOTO KAGAKU 🔿 KYOTOKAGAKU AMERICA INC. 63