Bathing and In-patient Care Simulator (Male)

Elderly whole-body manikin for basic patient care

This manikin facilitates training for various procedures essential for elderly patient care, such as bathing, hair care, positioning and oral care.



FEATURES

- 1 | Flexible and seamless joints
- 2 | Waterproof body
- 3 | Enables practices for shampooing, brushing and drying patient's hair
- 4 | Training for patient handling such as positioning, carrying or transferring on and off the bed





Q

Kyoto Kagaku Product Lineup Web



SPECIFICATIONS

SKILLS

Size: 168cm/ 66 in Weight: 14 kg / 30 lb Soft resin SET INCLUDES Latex free

Positioning

Hair care Oral care

Bathing

Patient handling

Range of joint-motion training

Diaper changing, clothing

1 manikin

1 denture

- 1 wig
- 1 underpants

30 lb Soft resin Latex free

MATERIALS

DESCRIPTIONS REPLACEMENT PARTS 11274-020 wig for M72