### M93C

# CVC Insertion Simulator III

# CVC III facilitates a sequence of uninterrupted simulation experience for specialist training

From aseptic field preparation to final catheter placement









#### FEATURES

- 1 Insertion practice is repeatable to obtain the best technique for patient safety.
- 2 |Both landmark guided insertion and ultrasound guided techniques can be practiced.
- 3 |The frictionless tissue-like simulation pad to facilitate the realistic sensation of CVC placement.

## SKILLS

Take complete steps to a successful CVC insertion

 Perform hand hygiene practices and prepare aseptic field
Decide on insertion site via ultrasound guidance or observation and palpation of landmarks

- | Insertion under ultrasound guidance or via landmark technique | Confirm flashback
- | Advance guidewire and sheath
- | Insert catheter





Pad material simulates features of human tissue



Anatomically accurate features including 3 insertion sites and rib cage

SPECIFICATIONS Size: W40 x D20 x H 34

#### Size: W40 x D20 x H 34 cm W15.7 x D8 x H13.5 inch SET INCLUDES

- 1 male upper torso manikin
- 2 CVC placement pad
- 1 skin for cannulation training
- 1 transparent anatomical block
- 2 artery tube (transparent): 2 type
- 2 vein tube (blue): 2 type2 vein pipe
- 1 pulsation unit (for landmark technique)
- 1 irrigation bottle (for UG technique)
- 1 plastic jar

# DESCRIPTIONS MATERIALS

Soft resin, hydrogel, hard resin Latex free

sample needle
syringe (2.5ml, 50ml each)

- 1 introductory ultrasound training block
  - REAL VESSELS

#### RECOMMENDED DEVICES Ultrasound machine and needles are r

Ultrasound machine and needles are not included.

#### REPLACEMENT PARTS

11347-320 irrigation bottle

11347-210 2 introductory ultrasound training block REAL VESSELS 11347-240 2 CVC placement pads 11347-260 1 pulsation unit (for M93UB)

🕻 KYOTO KAGAKU 🏷 KYOTOKAGAKU AMERICA INC. 58

Kyoto Kagaku Product Lineup Web