Pitting Edema Simulator

Five variations of pitting edema simulator including a normal one

Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"

Product Supervision

Junji Otaki Professor of Tokyo Medical University







FEATURES

- 1 | Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)"
- 2 |The models can be attached to the lower leg of simulated patients or training manikins.
- 3 4 levels of pitting edema according to pressures and the depth of indentations.
- 4 Repeatable and accurate, as if palpating actual skin surface to the fibula.

SKILLS

| Differentiation of the pitting edema:

- four different levels of edemas and a normal case
- | Assessment of edemas on the lower legs



SPECIFICATIONS

Size: 10x20cm / 3.9x7.9 inch Weight: 0.15kg/ 0.33lbs

SET INCLUDES

5 variations of strap-on pitting edema model 1 piece each

1 instruction manual

DESCRIPTIONS MATERIALS Soft resin Latex free

🕻 KYOTO KAGAKU 🥎 KYOTOKAGAKU AMERICA INC. 38

Q