MW24

3D Colonoscope Training Simulator NKS

Three dimensional representation of colon, based on analytical study of CT Colonoscopy (CTC) data

The three dimensional structure of the colon provides hands on training in a variety of skills such as prevention of loop forming, reduction of loops including a loop at LCF as well as insertion technique which forms a loop intentionally.

Product Supervision

Katuhisa Noda, MD, Noda-Genki Clinic Takatoshi Kitada, MD, Kitada Clinic Yasumoto Suzuki, MD, Matsushima Clinic





N loop

FEATURES

- 1 | Transparent skeleton body allows direct observation of the position of the tip of the scope as well as form of the colon.
- 2 | Stable rectal anatomy: the simulator is designed to acquire colonoscopy skills that use three Houston's valve as landmarks, which appear at the same position every time.
- 3 The sigmoid colon can be pre-set to take up any one of the three most common morphologies, which are deprived from analysis of CTC. The patterns can be changed quickly by sliding attachments and moving the colon tube.
- 4 The simulator is entirely water-resistant and can be cleaned easily by water.
- 5 | Light weight and compact, facilitating handling and transportation.

SKILLS

- Colonoscope insertion in three fundamental patterns: -Short alpha loop (pull back technique) -Long alpha loop (push forward technique)
- -N loop (pull back/push forward technique) | Colonoscope insertion in three patient positions:
 - -Left lateral
 - -Right lateral
 - -Spine
- | Manual abdominal compression

SPECIFICATIONS

Size: 43x28x16cm/17x11x6.3 inch Weight: 6kg/13.2lbs

MATERIALS

Soft resin, hard resin Latex free

SET INCLUDES 1 skeleton body

- abdominal membrane 1 body base
- colon-rectum tube
- (with a cap and a plug)

DESCRIPTIONS

- 1 slip resistant mat
 - plastic jar
 - 15 lubricant
 - 15 colonoscope gel
 - 1 storage case
 - 1 instruction manual

REPLACEMENT PARTS 11403-010

colon-rectum tube with a cap and a plug for MW24

11403-020

lubricant packs blue for MW24 (15x50ml)

- 11403-030 colonoscope gel for M24 (20ml each)

* Endoscope is not included.

🕻 KYOTO KAGAKU 🔿 KYOTOKAGAKU AMERICA INC. 55

Kyoto Kagaku Product Lineup Web