MW29

Oral Care Simulator Advanced



Oral care simulator which can mimic various patient conditions

True-to-life oral cavity and denture sets Actual water can be used to clean oral cavity



FEATURES

- 1 Viewable and accessible trachea and esophagus
- 2 Interchangeable denture with tooth diseases
- 3 Simple set-up and head positioning
- 4 | Includes washable plague and oral thrush
- 5 | Suction training with water







SKILLS

| Oral cavity assessment

- Care of the oral mucous membrane
- Brushing and cleaning of the oral cavity
- | Identifying oral thrush
- Care for humidity retention
- Replacement and maintenance of the denture
- Understanding of airway suction
- Massage of oral cavity

KEY FEATURES

Oral anatomy of denture with problems

spaced arch food residue anodontia stomatitis isolated tooth oral thrush

crowding gingival recession | root fracture | stump of a tooth



Patient conditions

Complete denture | Dental problems Toothless Intratracheal intubation | Tracheostomy | Tube feeding

Patient positions

Sitting position Half sitting position | Lateral position

SPECIFICATIONS

Size: W21xD37xH26 cm / 8.57x12.6x10.6inch Weight: 2.6kg/5.7lbs

SET INCLUDES

- 1 head/neck manikin
- base 1
- normal denture 1 denture with tooth disease
- 1 denture case
- sponge-tipped swabs toothbrush

2

mirror

1 simulated residue

- 1 denture brush
- 1 storage case
 - instruction manual

DESCRIPTIONS

MATERIALS Soft resin

Latex free

REPLACEMENT PARTS

- 11408-020 normal denture
 - 11408-030 denture with tooth disease
 - 11408-010 simulated residue
 - * Instruments are not included.

KYOTO KAGAKU 🔿 KYOTOKAGAKU AMERICA INC. 29

