Pressure Sores Simulator MW63 11446-000 Full Set 11446-100 **Basic Set** Set includes 11446-000 11446-100 Lower body torso Sacral pressure sore pad Stage II Stage III \bigcirc \bigcirc Stage IV (with granulation) Stage IV (with pockets) Sciatic pressure sore pad stage IV \bigcirc Two types of echo pads Main body Necrotic tissue removal set -Main body : Adult lower body torso Consumables set (simulated pus, blood, exudate and necrotic tissue) \bigcirc \bigcirc • Size : $W44 \times D40 \times H38$ Material : Soft resin, Latex free \bigcirc Storage case Weight: 3.8kg \bigcirc \bigcirc Instruction manual







simulated exudate simulated blood





Consumables set

simulated pus simulated exudate simulated blood simulated necrotic tissue



Necrotic tissue removal set



simulated pus





simulated necrotic tissue

Base for necrotic tissue removal Pads for necrotic tissue removal

Two types of echo pads (normal / edema)

11446-000

11446-100

KYOTO KAGAKU co.. LTD WORLDWIDE INQUIRIES & ORDERING

KYOTOKAGAKU HEAD OFFICE & FACTORIES

15 Kitanekoya-cho, Fushimi-ku, Kyoto 612-8388, JAPAN



http://www.kyotokagaku.com

w rw-kyoto@kyotokagaku.co.jp



http://www.kkamerica-inc.com/

🔀 info@kkamerica-inc.com

KYOTO KAGAKU

MW63 Pressure Sores Simulator



2021.08 G. 3000

OIN,

Assessment of Pressure Sores

A. Sacral pressure sore

This simulator includes a set of sacral pressure sore pads with different stages (4 types)





Partial-thickness loss of dermis Full-thickness skin loss



Full-thickness skin loss with granulation



B. Greater trochanter pressure sore C. Sciatic pressure sore



Full-thickness skin loss



Full thickness tissue loss with exposed bone

OIN,

Assessment of deep tissue by ultrasound

Sciatic pressure sore:

Ultrasound assessment of deep tissue injury Two types are included : normal and edema





edema



In accordance with the DESIGN-R 2020, all seven items can be trained using a combination of different pads and simulated liquid.

Depth

Exudate

Size

nflammation / nfection

Granulation

Neocrotic tissue

Pocket

*Japanese Society of Pressure Ulcers



0 III

Prevention and treatment of pressure ulcers / Cleansing /Dressing





cleansing

0 W

applying dressings

Realistic reproduction of human tissue allows for accurate training environment





	simulated blood	simulated exudate
Product Supervision	Masaharu Sugimoto, Part-	Terashi , Plastic and Reconstructi time Lecturer, Toyama University 1ent, Kobe University Hospital
SKILLS	Palpation of the wou	ent of deep tissue injury
	(edema in the sciatic Wound cleansing ca	

Applying dressings

NAL.N

Cleansing, applying dressings or topical medications and procedures of Negative Pressure Wound Therapy (NPWT)

Simulated exudate, and necrotic tissue set are included. Excellent quality of appearance and texture.

020

- Applying topical medications
- Negative pressure wound therapy
- Debridement of necrotic tissue
- Prevention of pressure sore