MW9

Peripheral Venous Catheter Placement Simulator

Holistic simulation of peripheral intravenous access management in just one arm



FEATURES

- 1 | IVC insertion and management from start to finish
- 2 | Two puncture sites: median antebrachial vein and dorsal veins
- 3 | Durable puncture pad

SKILLS

| Infection control Use of tourniquet | Identifying puncture site Sterilization of puncture site Puncture with IV cannula | Confirmation of flashback Application of pressure on vein Preparing and administering IV medications (simulated) | Dressing and securing cannula



Confirmation of flash back



Fixation



Application of pressure to vein

SPECIFICATIONS Size: W52 x D10 x H10 cm W20.5 x D3.9 x H3.9 inch

SET INCLUDES

- 1 arm model
- injection pads for dorsal veins of hand 2
- 2 injection pads for median antebrachial vein
- 1 bottle for simulated blood

- 1 stand for bottle of simulated blood
 - infusion set infusion bag
 - syringe (50ml)
- 1 plastic beaker

DESCRIPTIONS MATERIALS Soft resin

Latex free

10 simulated blood sticks (swab type) carrying bag 1

1 instruction manual

REPLACEMENT PARTS 11388-200 2 injection pads for median antebrachial vein 11388-300 2 injection pads for dorsal veins of hand

RECOMMENDED DEVICES

IV catheter needle 22G,

butterfly needle 23G

11388-400 10 simulated blood sticks (swab type)

KYOTO KAGAKU 🔿 KYOTOKAGAKU AMERICA INC. 66

1