C.H.A.R.L.I.E. Neonatal Resuscitation Simulator product profile

Product Name	C.H.A.R.L.I.E. Neonatal Resuscitation Simulator
Addresses Training Needs/Curricula	 Meets course curriculum for NRP[®] and Pediatric Advanced Life Support (PALS) Suitable for pre- and in-hospital training (EMS professionals and nurses)
Configuration	 Full-body manual simulator Multiple Packages for Optimal Training Heat-resistant (can be used in neonatal ICU unit and incubator) Features two birth anomalies for simulation SimVs offers an advanced simulation experience with C.H.A.R.L.I.E. NRP[®] Neonate with built-in scenario library and checklists
Dimension	 19 1/2" (49.53 cm) from crown to heel Approximately 8 1/2 lbs (3.86 kg)
Product Code/SKU	LF01420 CHARLIE with Interactive ECG LF01421 CHARLIE LF0142101 CHARLIE with Sim VS Pre-Hospital Defibrillator Simulation System LF0142102 CHARLIE with Sim VS Transport Simulation System LF0142103 CHARLIE with Sim VS Nursing Essentials Simulation System LF0142104 CHARLIE with Sim VS Nursing Med-Surg Simulation System



Nasco Healthcare 16 Simulaids Drive Saugerties, NY 12477 Phone: 1.833.NASCOHC (627-2642) Email: info@nascohealthcare.com www.nascohealthcare.com



Ready to save lives. Right from the start.

C.H.A.R.L.I.E. Neonatal Resuscitation Simulator:

C.H.A.R.L.I.E. Neonatal Resuscitation Simulator

Developed in conjunction with the latest Neonatal Resuscitation Program Guidelines (NRP[®])

- C: Compression and cardiac
- H: Heat compatible can be put in an incubator without melting
- A: Airway, arterial access via the umbilicus
- **R:** Resuscitation
- L: Laryngeal mask capable for oxygen use or PPV
- I: Intravenous, intraosseous, intubation
- E: ECG



Full-term newborn simulator for training on the most common neonatal care scenarios: Enables nurses and EMS workers to practice the essentials of neonatal care.



Add SimVs to upgrade to C.H.A.R.L.I.E. NRP[®] Neonate - an advanced simulation experience
Includes an instructor's tablet that controls the displays of 3 types of monitor interfaces: defibrillator, hospital monitor, and mechanical defibrillator

This manikin is intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.

