


PH-64 | 41928-000

## PET/ SPECT Brain Phantom



### Nucleus and the caudate nucleus for I-123 DatSCAN



#### FEATURES

- Anthropomorphic head phantom simulates the absorption and scatter characteristics of head and skull
- RI solution can be injected to striatum and putamen

#### APPLICATIONS

- SPECT
- DatSCAN

#### ANATOMY

Brain ventricle	Putamen
Striatum	Cerebrum
Caudate nucleus	Skull

	Volume
Cerebrum	1038.90 ml
Putamen (right)	5.46 ml
Putamen (left)	6.7 ml
Caudate nucleus (right)	5.38 ml
Caudate nucleus (left)	5.09 ml

#### DESCRIPTIONS

##### SET INCLUDES

- 1 head phantom (separated at the top)
- 1 brain container
- 1 storage case manual

##### MATERIALS

- Striatum case: epoxy resin
- Brain ventricle: polyurethane
- Brain case: polyurethane
- Artificial bone: epoxy resin
- Soft tissue: polyurethane

##### SPECIFICATIONS

- Phantom size:  
W20 x D21 x H33 cm  
W7.8 x D8.2 x H13 in
- Phantom weight:  
5.5 kg/ 12.1 lb

PH-53 | 41918-000

## Brain Phantom IB-20 advanced



### SHOW MORE!



#### APPLICATIONS

- SPECT, PET

EVALUATION PARAMETERS: common with IB-20

#### DESCRIPTIONS

##### SET INCLUDES

- 2 bone scatterer cases  
adult male: equivalent HU750  
elderly female: equivalent HU530
- 1 brain striatum phantom
- 1 screwdriver

- 1 velcro tape
- 1 petroleum jelly
- 1 storage case manual

##### SPECIFICATIONS

- Phantom size:  
W21 x D15 x H8 cm  
W8.2 x D5.9 x H3.1 in

- Brain striatum: epoxy resin
- Brain striatum container: urethane resin
- Cerebral ventricle: urethane resin
- Brain stratum phantom cover: acrylic
- Bone scatterer case: epoxy resin

PH-27 | 41530-000

## Brain Phantom IB-10



### SHOW MORE!



#### APPLICATIONS

- SPECT, PET

#### EVALUATION PARAMETERS

- Homogeneity evaluation
- Cross calibration
- Gamma ray absorption rate by a skull

- Detectivity of gray matter and white matter
- Spatial resolution of negative images (IB-10 set only)
- Radioactive concentration and linearity of SPECT value (IB-10 set only)

#### DESCRIPTIONS

##### SET INCLUDES

- 1 brain unit
- 1 skull container unit
- 1 J-Jack phantom
- 1 section phantom

##### MATERIALS

- Acrylic resin/ urethane resin

##### SPECIFICATIONS

- Phantom size:  
W21 x D15 x H8 cm/ W8.2 x D5.9 x H3.1 in