US-19

Lung Ultrasound Training Phantom ECHOLUNG

Simple, durable and effective tool to learn the basics of lung ultrasound.

Product Supervision Hirosaki University School of Medicine & Hospital Department of General Medicine Tadashi Kobayashi



FEATURES

- 1 | Four kinds of units: A-LINE, B-LINE(few), B-LINE (many), pneumonia and pleural effusion to learn the basics of lung ultrasound.
- 2 | Anatomy includes: lung, chest wall, ribs
- 3 No special storage needed
- 4 | Suitable for training in large groups

SKILLS

| Handling and manipulation of probe on lung ultrasound. | Scanning the lungs to evaluate pneumonia findings.



A-LINE



B-LINE



B-LINE(many)



Pneumonia and pleural effusion

SPECIFICATIONS

- Size:Body:32×25×12.5m/13×10×4.9 inch Module:14.8×11.8×11m/5.8×4.6×4.3 inch
- Weight:Body:1.9kg A-LINE/B-LINE/B-LINE (many): 0.9kg Pneumonia and pleural effusion: 1.4kg

DESCRIPTIONS -

SET INCLUDES

- 1 Main Body
- 4 Kinds of lung modules
- 1 Instruction manual
- 1 Talcum powder
- 1 Storage case

Q