

US-8




Pediatric FAST/Acute Abdomen Phantom

Ultrasound training for pediatric FAST procedures

This world's first pediatric ultrasound torso phantom provides opportunities of hands-on training in ultrasound examination of a two-year-old patient

Product Supervision

Shunsuke Nosaka, MD, PhD
Director of Radiology
National Medical Center for Children and Mothers




FEATURES

- 1 | Thoracoabdominal organs, bone structure and free-fluid for practicing FAST examination procedures
- 2 | Commonly seen pathologies in pediatrics
- 3 | Acquire fundamental pediatric abdominal ultrasound skills, as well as FAST/Acute sonography skills

SKILLS

- | FAST procedures:
Internal hemorrhage at perihepatic, perisplenic, pelvis, and pericardium area
- | Sonography for acute pediatric patient:
Bowel intussusception, appendicitis and biliary dilatation

KEY FEATURES



Pleural hemorrhage



Hydronephrosis



Appendicitis

DESCRIPTIONS

SPECIFICATIONS

Size: W41 x D17 x H15 cm Weight: 7kg / 15.4lbs
W16.1 x D6.7 x H5.9 inch

SET INCLUDES

- 1 | ultrasound phantom
- 1 | talcum powder
- 1 | storage case
- 1 | tutorial manual (DVD)

MATERIALS

Polyurethane elastomer
Latex free

REPLACEMENT PARTS

RECOMMENDED DEVICES

* Ultrasound scanner is not included.

