

M95

Simulator Pediatric IV Hand



Set of pediatric hands for practicing IV
Become confident in safe insertion of intravenous cannulas on small infant hands with thin and soft veins

FEATURES

- 1 | Fingers and dorsum of hand are flexible, able to be gripped, and repositionable as in real-life procedures.
- 2 | Two puncture sites on the back of the hand
- 3 | Practice dressing and securing the cannula

SKILLS

- | Pediatric dorsal veins injection
- | Blood collection
- | Intravenous infusion

DESCRIPTIONS

SPECIFICATIONS

Size: W18 x D30 x H9.5 cm / W7.1 x D11.8x 3.7 inch (one-year-old)
W18 x D36 x H10 cm / W7.1 x D14.2x 3.9inch (three-year-old)

Power: AC100-240V
Power Consumption: 5W

SET INCLUDES

2 injection arm-hands (one-year-old)	1 spoon
2 injection arm-hands (three-year-old)	1 plastic jar
2 stands	1 plastic beaker
1 circulation pump	1 storage case
1 cotton swab (red)	1 instruction manual

MATERIALS

Soft resin, hard resin
Latex free

REPLACEMENT PARTS

11350-010	1 replacement injection arm-hand (three-year-old)
11350-020	1 replacement injection arm-hand (one-year-old)
11388-400	10 cotton swab for coloring (red)

RECOMMENDED DEVICES

indwelling needle: 23G
butterfly needle: 23G

* Needles are not included.

M43D


Pediatric Lumbar Puncture Simulator II

Acquire skills for safe pediatric lumbar puncture

The full body is soft and bendable to facilitate realistic body positioning required for a safe and successful lumbar puncture.



Level of posterior superior iliac crests



Product Supervision

Takahiro Amano, M.D.
Vice President
Senior Vice-Dean, Postgraduate School
Professor and Director, Center of Postgraduate
Medical Education International University of Health
and Welfare,
Honorary Director, Sanno Medical Center
Gregory A. Plotnikoff, M.D., M.T.S., F.A.C.P.
Medical Director, Abbott Northwestern Hospital

FEATURES

- 1 | Soft, seamless baby manikin allows practice in safe patient positioning and holding.
- 2 | Two types of puncture pads: deep or shallow
- 3 | Realistic needle tip resistance
- 4 | Anatomically correct landmarks for palpation

SKILLS

- | Patient positioning
- | Palpation of landmarks
- | Lumbar puncture and CSF collection
- | CSF pressure measurement

DESCRIPTIONS

SPECIFICATIONS

Size: W43 x D34 x H20 cm
W17 x D13.4 x H7.8 inch

Weight: 1.5 kg / 3.3 lbs

SET INCLUDES

1 baby manikin (7-10 months)	1 syringe (50mL)
2 pediatric LP block (shallow)	1 instruction manual
2 pediatric LP block (deep)	
1 irrigation bag	
1 stand for irrigation bag	

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11348-100-04	stand for irrigation bag
11348-100-05	irrigation bag
11348-410	pediatric LP block (shallow)
11348-420	pediatric LP block (deep)

RECOMMENDED DEVICES

22G-24G spinal needle

* Needles are not included.