MW14

Difficult Airway Management Simulator -Demonstration Model-

Great addition to a demonstration tool for airway management

Compact table-top design, anatomical accurate and a variety of DAM settings. Ideal for demonstration of skills and devices.

FEATURES



Difficult Airway Management Simulator for Bronchofiberscopy

Realistic trainings in insertion and management of a fiberscope

Anatomically correct trachea and bronchi provides realistic view through a scope, allowing recognition of bifurcations.





FEATURES

1 | Anatomically correct trachea and bronchi that include trachea bifurcation for segmental bronchi.

SKILLS

Insertion and management of a fiberscope

- Recognition of bifurcations by a fiberscope
- Airway opening techniques (head tilt, jaw thrust)
- Sniffing position

SPECIFICATIONS

Size: W45 x D20 x H25 cm W17.7 x 7.9 x H9.8 inch Weight: 2.65kg/ 5.84lbs

SET INCLUDES

1 manikin

- 3 upper incisors
- 1 lubricant
- 1 syringe
- 1 carrying bag

1 instruction manual

ringe rrying bag

MATERIALS Soft resin Latex free

DESCRIPTIONS

REPLACEMENT PARTS

- 11393-020 1 face mask II 11392-080 2 replacement tongues II 11392-010 10 upper incisors
- 11229-050 lubricant

RECOMMENDED DEVICES

(Macintosh laryngoscope) blade size 3, 4 (Endotracheal tube) internal diameter 7.0mm, 7.5mm (Laryngeal mask) Air-Q: size 3.5 LMA Supreme: size 4 I-gel: size 4

* Instruments are not included.

KYOTO KAGAKU

99