

# PROMPT FLEX RANGE

Designed with, and used by leading worldwide Clinicians. The cost effective, modular nature of the **Birthing Simulator - PROMPT Flex** allows the training and practice of birthing, both normal and difficult, along with several maternal emergencies.

Used with Simulated patients, the key skill of communication can be practiced.

Multiple studies show positive outcomes through improved teamwork after a system wide commitment to enduring training.





Learning has never been more lifelike limbsandthings.com

# **STANDARD & ADVANCED BIRTHING SIMULATOR - PROMPT FLEX** ARE STAND-ALONE TASK TRAINERS THAT CONTAIN A BIRTHING MOTHER & BABY WITH KEY ANATOMICAL FEATURES











Standard only

### Birthing Simulator - PROMPT Flex -Standard

#### 80200

A modular design birthing simulator allowing for numerous training scenarios.

Used either for stand-alone training or integrated into hybrid simulation, it is an ideal training solution for all skills relating to routine and difficult deliveries.

## Birthing Simulator - PROMPT Flex -Advanced



Incorporates wireless force monitoring where force applied to the baby's head is measured during a shoulder dystocia drill. Includes an app with inbuilt scenarios allowing monitoring and recording of force, interventions made and time taken to deliver the baby.



Tablet not included Download app at limbsandthings.com





	80100	80106
SKILLS	80200	80206
• Communication & team work skills (SP)	$\checkmark$	~
• Normal	$\checkmark$	$\checkmark$
• Breech	$\checkmark$	$\checkmark$
<ul> <li>Shoulder dystocia management</li> </ul>	$\checkmark$	$\checkmark$
• Instrumental (forceps & vacuum devices	) 🗸	$\checkmark$
<ul> <li>Delivery of placenta</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Cord prolapse</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Urinary catheter placement</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>IM injection pads</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Force feedback of baby's head</li> </ul>	×	$\checkmark$
• App with scenarios & force monitoring	×	$\checkmark$

WHAT'S INCLUDED	80200	80206
Birthing Mother (+Upper Legs)		~
Standard Baby	$\checkmark$	×
Wireless Force Monitoring Baby	×	$\checkmark$
Placenta	$\checkmark$	$\checkmark$
Abdomen for PROMPT Flex	$\checkmark$	$\checkmark$
Perineum and Birth Canal	$\checkmark$	$\checkmark$
Bed Straps	$\checkmark$	$\checkmark$
Lubricant	$\checkmark$	$\checkmark$
Wheeled Carry Case	$\checkmark$	$\checkmark$

80100 80106

#### QUALITIES

#### ANATOMY

- Birth canal and cervix
- Ischial spines and pubic bone
- Gynaecoid pelvis
- Articulating thighs
- IM injecting
- Fully articulated baby with detachable umbilical cord and placenta
- Baby includes clavicles, fontanelles and flexible head

#### **KEY FEATURES**

- Improved and more realistic pelvic floor
- Articulating thighs for McRoberts procedure
- Stretchable perineum
- Soft flexible birthing canal
- Additional modules to extend learning
- Optional lower legs for extended realism
- Optional wireless force monitoring for training in shoulder dystocia management

#### SAFETY

• Latex free

#### **CLEANING & CARE**

• Skin is washable using soap and water

# Cervical Cerclage Module (CCM) - PROMPT Flex

#### 80280

This module is the ideal tool for trainees looking to develop the skills of elective and emergency cervical cerclage.

The easy to use design offers realistic movement of the cervix down the vaginal vault, enabling trainees to place the cerclage as high up the cervix as possible, which supports the learning of the McDonald technique.



Elective

#### Emergency

#### SKILLS

- · Management of cervical incompetence
- Elective cerclage
- Emergency (or rescue) cerclage
- McDonald cerclage technique

#### QUALITIES

#### REALISM

- Realistic representation of a dilating cervix with presenting amniotic sac for emergency procedures
- Restrictive space within the vagina to perform the procedure as in life

#### **KEY FEATURES**

- Detachable cervix to allow for repeated training in both elective and emergency cerclage procedures
- The detachable OS provides instant feedback to the trainee
- Fluid filled amniotic sac replicates bulging membranes for emergency procedures

#### VERSATILITY

- Allows for the practice of both elective and emergency cervical cerclage
- Can be quickly and easily fitted into PROMPT Flex

#### SAFETY

• Main product is latex free however Amniotic Sacs provided contain latex

Developed in collaboration with:

- Dr Graham Tydeman, Consultant Obstetrician and Inventor, NHS Fife
- Professor Andrew Shennan OBE and Professor Annette Briley, Guy's and St. Thomas' NHS Foundation Trust



# WHAT'S INCLUDED

Cervical Cerclage Frame Cervical Cerclage Vagina

- Cervical Cerclage Elective OS (x5)
- Cervical Cerclage Emergency OS (x5)
- Cervical Cerclage Amniotic Sacs (x5)
- Cervical Cerclage Fluid System
- Cervical Cerclage Perineum

# Cervical Dilatation & Effacement Module (CD&E) - PROMPT Flex

80102 80202

Allows for training in assessment of both the latent and active first stages of labour.

This versatile module is an optional add-on to the PROMPT Flex Standard or Advanced.







#### REALISM

• Realistic representation of cervixes, including anterior lip, and presenting parts in soft birth canal, with palpable ischial spines

#### **KEY FEATURES**

- Inserts to represent early labour cervixes effacement, dilation and ripeness in line with Bishop's scoring
- Positioning mechanism allows adjustment of station, dilation and tilt without removal from the simulator
- Dynamic positioning mechanism allows adjustment of dilation in active labour
- Markers allow tutor to read positioning in situ
- Numerous presenting part inserts including
  - flexed
  - deflexed
  - brow
  - face
  - breech
  - caput and moulding

#### VERSATILITY

• Presenting parts and cervix can be set up outside the model for demonstration, then placed inside for examination

#### SAFETY

• Latex free











Fetal station adjustment (-3 to +3)

Cervical Dilation (1-10cm)

### SKILLS

- Assessment and Bishop's scoring of:
  - Cervical dilation (1-10cm)
  - Cervical effacement (0-100%)
  - Cervical ripeness/consistency (soft, medium, hard)
  - Cervical position (anterior, mid, posterior)
  - Fetal station (-3 to +3)
- Assessment of and artificial rupture of membranes
- Assessment of presenting part flexed, deflexed, brow, face, breech, caput and moulding and caput formation

#### WHAT'S INCLUDED

#### 0102 80202

Perineum
Pelvic Ring
Static Rig
Dynamic Rig
Set of 7 Latent Labour Cervixes
Dynamic Cervix
Set of 5 Presenting Parts
Set of 2 Difficult Labour Heads
Amniotic Membrane Set
Set of 3 Caputs (x2)

# Enhanced Caesarean Section Module (ECSM) - PROMPT Flex

80300 80400

Ideal for practising the necessary skills required when performing both routine and complex caesarean section. Trainees can have exposure to delivering a difficult baby such as an impacted fetal head (IFH), or a transverse lie.

A pre-incised abdomen allows for repeated delivery practice to help improve confidence and familiarity of manoeuvres.



#### QUALITIES

#### REALISM

- Abdominal skin has realistic appearance, fat layer and feel of at-term abdomen
- Realistic feel of uterus interior during delivery
- Head and shoulders present at realistic level in relation to the uterine incision

#### **KEY FEATURES**

- Uterus holds baby firmly in position for cephalic, breech and transverse presentations
- Fundal pressure can be applied to aid delivery of the baby for IFH manoeuvres
- Pre-incised skin allows easy repeat delivery practice
- Realistically mimics the restrictive space of the pelvic inlet for IFH manoeuvres
- Accommodates retraction using metal, plastic or rubber retractors

#### VERSATILITY

• Allows for training in routine and complex caesarean section at full dilation

#### SAFETY

• Latex free

The module can be used to practise any of the following methods of delivering a deeply impacted fetal head:

- Manual cephalic extraction manoeuvres
- Vaginal push up
- Reverse breech extraction (pull technique)
- Patwardhan's technique
- Use of a fetal pillow



#### SKILLS

- Cephalic and breech delivery
- C-section at full dilation:
  - Head disimpaction using vaginal push methods
  - Use of balloon disimpacting device
  - Reverse breech extraction
  - Patwardhan's techniques
- Transverse lie
- Instrument assisted C-section delivery
- Delivery of the placenta for routine c-sections
- Teamwork and human factors training in emergency C-section

#### WHAT'S INCLUDED

#### 80400

Perineum

- Uterus for Caesarean Section Module
  - Impacted Fetal Head Uterus
- Pre-incised Delivery Skin
  - Cervix for Caesarean Section Module
- Foam Uterus Support
  - Pelvic Floor

## Postpartum Haemorrhage Module - PROMPT Flex

80101 80201

This module is an optional add-on for the Standard or Advanced Birthing Simulators.

It allows training in the management of post-partum bleeding, and offers a realistic simulation of an atonic and contracting uterus.





#### QUALITIES

#### REALISM

• Atonic uterus can simulate contraction on fundal massage

#### **KEY FEATURES**

- Blood loss of up to 2 litres can be simulated
- Air bulbs are used to stimulate blood flow and control atonic state of the uterus

#### VERSATILITY

• Can be quickly and easily fitted into the birthing mother to allow simulation of situations where bleeding occurs directly after delivery

#### SKILLS

- Management of postpartum bleeding:
  - Fundal massage
  - Bi-manual compression
  - Haemostatic balloon insertion and management
- Communication skills with mother when used with Simulated Patient
- Intramuscular injection
- Estimating blood loss



#### WHAT'S INCLUDED

80201

PPH Uterus Assembly - PROMPT Flex Concentrated Venous Blood - Starter Pack (x2) PPH Blood Container - PROMPT Flex Using PROMPT in conjunction with effective professional training, reports show:

- 50% Reduction on brain injury
- 45% Reduction in school-aged cerebral palsy
- Reduction in maternal deaths

Source: <u>https://www.promptmaternity.org/publications</u>

#### R **PROMPT Flex** is suitable for use on Obstetric Training Programmes including: ACOG Course – ECO (Emergencies in Clinical Obstetrics) MORE OB (Managing Obstetric Risks Efficiently) ALARM (Advances in Labor and Risk Management) - Canada PROMPT (PRactical Obstetric Multi Professional Training) ALSO (Advanced Life Support in Obstetrics) **PRONTO International** BEOC (Basic Essential Obstetric Care) **REOT (Rural Emergency Obstetrics Training)** BLSO (Basic Life Support in Obstetrics) ROBUST (RCOG Operative Birthing Simulation Training) SCOTTIE (Scottish Core Obstetric Teaching and Training In Emergencies) Care Team OB (American Academy of Family Physicians) TeamSTEPPS (Team Strategies and Tools to Enhance Performance CEOC (Comprehensive Essential Obstetric Care) and Patient Safety) MOET (Managing Obstetric Emergencies and Trauma)

# YOU MAY ALSO BE INTERESTED IN CLEAN BLEED<sup>TM</sup> MAT

#### 70400

Ideal product to support any simulation scenario where external blood or fluid loss is key indicator and measure of emergency interventions taken.



### Designed to support different scenarios



TASK TRAINER







REQUEST A QUOTE



REQUEST A DEMO



INT V.10 Jan 2022

Limbs & Things

E: sales@limbsandthings.com T: 00 44 (0)117 311 0500

# limbsandthings.com