

MW48



Full Perinatal Patient Care Simulator “Hana”

Be a part of Hana's journey as a new mother


Training manikin from gestational to post-natal period
Let us help you take your first step in perinatal care!



FEATURES

- 1 | Interchangeable abdomen to simulate perinatal examination, ultrasound examination of 23-week-fetus, cervical examination, labor, perineal tear and puerperal.

POINTS




One-touch interchangeable parts



Compatible with birthing chairs




Durable yet realistically elastic



Accurate anatomy

SKILLS



Perinatal examination

During pregnancy

Perinatal examination

- Palpation
(e.g. using Leopold's maneuvers, cephalopelvic disproportion assessment)
- Repositionable fetus
- Auscultation of fetal sounds, with adjustable heart rates from 60 to 180
- Measurable abdomen girth, fundal height and external pelvimetry



Ultrasound examination

- Screening of 23-week-fetus
- Fetal size assessment using BPD, AC, FL, and so on
- Measurement of amniotic fluid volume
- Determination of fetus presentation



Childbirth

Cervical examination

- 5 types of cervical dilation models
- Determination of bishop score (cervix score)
- Check fetal head position




During labor

- Various positions: Supine, lateral, all-fours
- Perineal protection
- Cephalic and breech presentations
- Obstetric forceps delivery and vacuum extraction delivery
- Tying and cutting of the umbilical cord
- Delivery of the placenta
- Inspection of the amniotic sac
- Insertion of urinary catheter



Perineal tear

- Perineorrhaphy:
 - single interrupted suture
 - vertical mattress suture
- Repetitive suture practice for first degree perineal tear



Puerperal

Puerperal assessment

- Assessment uterine involution in early puerperium: 4 types of variations
- Measurement and palpation of uterine fundus
- Perineal cleaning and assessment
- Postpartum massage
- Assessment of anal prolapse



DESCRIPTIONS

SPECIFICATIONS

Size: 164 cm/ 63 inch Weight: 20 kg / 2.2 lbs

SET INCLUDES

* Components for basic manikin. Training modules are to be added according to requirements

1 manikin	1 injection pad (median vein)
1 wig	1 irrigation bag
1 face mask for Hana	1 drainage pump /drainage hose
1 denture	1 lubricant for training models
1 female genitalia	1 talcum powder
1 maternity clothes	instruction manual

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11229-050	lubricant for training models
11251-030	2 shoulder IM pads
11416-040	5 velamen sheets
11416-030	10 umbilical cords for omphalotomy
11417-010	2 genitalia insert
11415-040	lubricant (500ml)

RECOMMENDED DEVICES

-	
11406-040	wig
11415-030	genitalia unit
11415-010	cervical dilation inserts
11416-020	placenta